



Innomed eNews – Februar 2024

Basic Screw Removal System

System designed to help remove damaged and broken screws from 1.5 to 7.0 mm

Screw Removal Pliers



#2022-01

Mini Lexer Gouges



4 mm Gouge #2022-02



6 mm Gouge #2022-03



10 mm Gouge #2022-04

Sharp Hook



#2022-SH

T-Handle with AO-End



#2022-T

Extraction Screws



For 1.5/2.0 mm #2022-05



For 2.7/3.5/4.0 mm #2022-06

Set in Case



Product No's:

2022-00 [Complete Set]

Set Includes:

2022-01 [Screw Removal Pliers]

Overall Length: 8" (20,3 cm)

 For 4.5/5.0/6.5/7.0 mm #2022-07

2022-02 [Mini Lexer Gouge – 4 mm]
Overall Length: 7" (17,8 cm)
Gouge Width: 4 mm

Extraction Bolts

 For 1.5 mm Screw #2023-01

2022-03 [Mini Lexer Gouge – 6 mm]
Overall Length: 7" (17,8 cm)
Gouge Width: 6 mm

 For 2.0 mm Screw #2023-02

2022-04 [Mini Lexer Gouge – 10 mm]
Overall Length: 7" (17,8 cm)
Gouge Width: 10 mm

 For 2.7 mm Screw #2023-03

2022-05 [Extraction Screw for 1.5/2.0 mm]
Overall Length: 1.6" (4 cm)

 For 3.5/4.0 mm Screw #2023-04

2022-06 [Extraction Screw for 2.7/3.5/4.0 mm]
Overall Length: 1.6" (4 cm)

 For 4.5 mm Screw #2023-05

2022-07 [Extraction Screw for 4.5/5.0/6.5/7.0 mm]
Overall Length: 1.6" (4 cm)

 For 5.0/6.5/7.0 mm Screw #2023-06

2022-CASE [Case for Basic Screw Removal System]

Trephines

2022-IP [Instruction Plate]
Dimensions: 7.875" x 3.65" (20 x 9,3 cm)

 For 1.5 mm Screw #2023-07

2022-SH [Sharp Hook]
Overall Length: 6.1" (15,5 cm)

 For 2.0 mm Screw #2023-08

2022-T [T-Handle with A0-End]
Overall Length: 5.9" (15 cm)
Handle Width: 3.15" (8 cm)

 For 2.7 mm Screw #2023-09

 For 3.5/4.0 mm Screw #2023-10

2023-01 [Extraction Bolt for 1.5 mm Screw]
Overall Length: 2.35" (6 cm)

 For 4.5 mm Screw #2023-11

2023-02 [Extraction Bolt for 2.0 mm Screw]
Overall Length: 2.35" (6 cm)

 For 5.0/6.5/7.0 mm Screw #2023-12

2023-03 [Extraction Bolt for 2.7 mm Screw]
Overall Length: 2.35" (6 cm)

Spare Trephine Cutting Ends

2023-04 [Extraction Bolt for 3.5/4.0 mm Screw]
Overall Length: 2.35" (6 cm)

 For 1.5 mm Screw #2024-01

2023-05 [Extraction Bolt for 4.5 mm Screw]
Overall Length: 3.15" (8 cm)

 For 2.0 mm Screw #2024-02

2023-06 [Extraction Bolt for 5.0/6.5/7.0 mm Screw]
Overall Length: 3.94" (10 cm)

 For 2.7 mm Screw #2024-03

2023-07 [Trephine for 1.5 mm Screw]
Overall Length: 4.125" (10,5 cm)

 For 3.5/4.0 mm Screw #2024-04

2023-08 [Trephine for 2.0 mm Screw]
Overall Length: 4.125" (10,5 cm)


 For 4.5 mm Screw #2024-05

2023-09 [Trephine for 2.7 mm Screw]
Overall Length: 4.125" (10,5 cm)

 For 5.0/6.5/7.0 mm Screw #2024-06

2023-10 [Trephine for 3.5-4.0 mm Screw]
Overall Length: 4.125" (10,5 cm)

Instruction Plate #2022-IP



2023-11 [Trephine for 4.5 mm Screw]
Overall Length: 5.4" (13,7 cm)

2023-12 [Trephine for 5.0-7.0 mm Screw]
Overall Length: 5.4" (13,7 cm)

2024-01 [Spare Trephine Cutting End for 1.5 mm Screw]
Overall Length: 1.6" (4 cm)

2024-02 [Spare Trephine Cutting End for 2.0 mm Screw]
Overall Length: 1.6" (4 cm)

2024-03 [Spare Trephine Cutting End for 2.7 mm Screw]
Overall Length: 1.6" (4 cm)

2024-04 [Spare Trephine Cutting End for 3.5/4.0 mm Screw]
Overall Length: 1.6" (4 cm)

2024-05 [Spare Trephine Cutting End for 4.5 mm Screw]
Overall Length: 2.75" (7 cm)




2024-06 [Spare Trephine Cutting End for 5.0/6.5/7.0 mm Screw]
Overall Length: 2.75" (7 cm)

Carpal Tunnel Release Guide and Blade Set

Guide designed to help protect the median nerve while providing a track that allows for the smooth advance of the blade to divide the transverse carpal ligament during a mini-open, non-endoscopic approach

Guide designed by Peter J. Evans, MD, PhD





Product No's:

1124-00 [Set of One Guide and One Blade]

Set Includes:

[Carpal Tunnel Release Blade]
Overall Length: 4.75" (12 cm)
Blade Length: 1.375" (3,5 cm)
Blade Height: .22" (5,6 mm)
Blade Width: 1 mm

1128 [Evans Carpal Tunnel Guide]
Overall Length: 8" (20,3 cm)
Blade Guide Widths: 2 mm and 5 mm


Replacement Parts:

1124-02 [Carpal Tunnel Release Blades] **Pack of 2**
One blade included in set, two with this number

1128 [Evans Carpal Tunnel Guide]

Cobra Retractors with Blunt Teeth

General purpose hip instruments for use around the femur and acetabulum with teeth to help prevent slippage



Product No's:

6130-01 [Standard with Blunt Teeth]
Overall Length: 12" (30,5 cm)
Handle Length: 7" (17,8 cm)
Blade at Widest: 33 mm

6130-02 [Straight Tip with Blunt Teeth]
Overall Length: 12" (30,5 cm)
Handle Length: 7" (17,8 cm)
Blade at Widest: 33 mm

Chandran Thigh Lift Positioner

Designed to help lift and position the thigh from above during knee



Product No's:

4167-00 [Positioner Set]

Included in set and also sold Individually:

4167-01 [Cradle]

Dimensions: 7.5" x 7.75" (19,1 x 19,7 cm)

4167-02 [Rod]

Width: 13.6" (34,5 cm)

Height: 20.6" (52,3 cm)

4150-PS [Post Screw]

Two included in set, one with this product number

2760-P [Silicone Pad]

Dimensions: 9" x 8" x 1/2" thick (22,9 x 20,3 x 1,3 cm)

[Orange Strap]

Dimensions: 2" x 31.5" (5,1 x 80 cm)

Two included in set, see below for replacement 10-pack

Replacement Items:

8130-P [Orange Straps - 10-Pack]

Dimensions (Each): 2" x 31.5" (5,1 x 80 cm)

2760-P [Silicone Pad]

Dimensions: 9" x 8" x 1/2" thick (22,9 x 20,3 x 1,3 cm)

Optional Items (Not included in Set):

4167-03 [Adapter]

Width: 12.4" (31,5 cm)

Height: 4.875" (12,4 cm)

9125 [Rotating Table Clamp]

The optional thigh lift adapter is designed for use with a hydraulic lift device instead of the manual lift rod with table clamp.



Für weitere Fragen zu unseren Produkten stehen wir Ihnen gerne zur Verfügung!

Do not hesitate to contact us for further information about our products!

Tel. (for EU states): +49 (0) 7720 46110 60

Tel. (outside the EU): +41 (0) 41 740 67 74

For orders/offers: orders@innomed-europe.com

For loaners: loaners@innomed-europe.com

General Information: info@innomed-europe.com

[VISIT OUR WEBSITE](#)

Copyright © 2023 Innomed-Europe LLC, All rights reserved.

Sie erhalten diesen Newsletter da Sie als Kunde bei Innomed-Europe registriert sind.

Our mailing address is:
Innomed-Europe LLC
Alte Steinhauserstrasse 19
6330 Cham Switzerland

Want to change how you receive these emails?
You can [update your preferences](#) or [unsubscribe from this list](#).