




## Innomed eNews – März 2024

### Goytia Stackable Hohmann Retractors

*Interlocking design helps to increase depth and leverage in hip exposure, particularly of the anterior acetabulum—especially useful with large patients*



Designed by Robin N. Goytia, MD



Custom fitted holes for interlocking retractors helps provide stability.



When "stacked", the increased lever arm of the retractor helps reduce fatigue.

Ideal for use with large patients where extra depth, leverage and force is needed.

New deep versions are 2" (5 cm) deeper for use with large patients where extra depth, leverage and force is needed.

Product No's:

4551-D [Deep Standard]  
Overall Length: 11.4" (29 cm)  
Handle Length: 5.625" (14,3 cm)  
Blade Width: 20 mm

4552-D [Deep Bent]  
Overall Length: 9.1" (23,1 cm)  
Handle Length: 5.625" (14,3 cm)  
Blade Width: 20 mm

4553-D [Deep Wide]  
Overall Length: 12" (30,5 cm)  
Handle Length: 5.625" (14,3 cm)

Blade Width: 4.3 mm

### Johnson Low Profile Foot & Ankle Retractors

*Designed for soft tissue retraction in the foot and ankle*



Designed by Michael Johnson, MD



Product No's:

1636-01 [Straight]  
Overall Length: 7.5" (19 cm)  
Blade Width: .41" (1 cm)

1636-02 [Bent]  
Overall Length: 8.875" (22,5 cm)  
Blade Width: .41" (1 cm)

### Mueller-type Femoral Neck Elevator with Blunt Teeth

*Designed to elevate the proximal femur, the additional blunt teeth on the end allow for better gripping*



Product No's:

3415-01  
Overall Length: 13.5" (34,3 cm)  
Handle Length: 5.75" (14,6 cm)

### Proximal Tibia Resection Shark Tooth Clamp

*Designed to help grasp and remove the cut proximal portion of the tibia during total and uni knee arthroplasty*



Designed by Shara Diers, PA-C



Product No: 3651  
Overall Length: 7" (17,7 cm)  
Platform Width: .79" (2 cm)

Für weitere Fragen zu unseren Produkten stehen wir Ihnen gerne zur Verfügung!  
Do not hesitate to contact us for further information about our products!

**Tel. (for EU states): +49 (0) 7720 46110 60**

**Tel. (outside the EU): +41 (0) 41 740 67 74**

**For orders/offers: [orders@innomed-europe.com](mailto:orders@innomed-europe.com)**

**For loaners: [loaners@innomed-europe.com](mailto:loaners@innomed-europe.com)**

**General Information: [info@innomed-europe.com](mailto:info@innomed-europe.com)**

**VISIT OUR WEBSITE**

*Copyright © 2023 Innomed-Europe LLC, All rights reserved.*

Sie erhalten diesen Newsletter da Sie als Kunde bei Innomed-Europe registriert sind.

Our mailing address is:  
Innomed-Europe LLC  
Alte Steinhauserstrasse 19  
6330 Cham Switzerland

Want to change how you receive these emails?  
You can [update your preferences](#) or [unsubscribe from this list](#).