




Innomed eNews – Juli 2024

Offset Canal Finder

A smooth, double bent, offset canal finder designed to assist during anterior hip surgery



Product No: 4987
Overall Length: 14" (35,6 cm)
Handle Length: 5.375" (13,6 cm)

Offset Gouge for Posterior Osteophyte Removal in TKA

An offset gouge with a rounded edge designed to more effectively remove osteophytes from round posterior condyles during total knee arthroplasty



Designed by Robert Steensen, MD

Product No: 3731
Overall Length: 11" (27,9 cm)
Blade Width: .815" (2 cm)



PCL Retractor – Straight

Designed to straddle the cruciate ligament



Designed to straddle the cruciate ligament and lie in the femoral condylar notch, allowing the surgeon to retract the tibia away from the femur for better access.

Product No: 2820-S [Straight]
Overall Length: 10.375" (26,3 cm)
Prong Width: 5 mm | 10 mm Gap | 5 mm

Powers Flared Small Knee Retractor

A bent knee retractor with a cobra flare to help provide optimal exposure



Designed by Mark Powers, MD

Product No: 6291
Overall Length: 8" (20,3 cm)
Pad Flare Width: 17 mm
Point Length Below Pad: 21 mm

Rogzinski Glenoid Reaming Retractor



Product No: 4277-01
Overall Length: 11.25" (28,6 cm)
Handle Width: 4.5" (11,4 cm)



Screw Removal Locking Pliers

Can be used to help remove bone from around screw heads or broken screws



Product No: 2022-01
Overall Length: 8" (20,3 cm)

T-Handle with A0-End



Product No: 2022-T
Overall Length: 5.9" (15 cm)
Handle Width: 3.15" (8 cm)

Für weitere Fragen zu unseren Produkten stehen wir Ihnen gerne zur Verfügung!
Do not hesitate to contact us for further information about our products!

Tel. (for EU states): +49 (0) 7720 46110 60

Tel. (outside the EU): +41 (0) 41 740 67 74

For orders/offers: orders@innomed-europe.com

For loaners: loaners@innomed-europe.com

General Information: info@innomed-europe.com

VISIT OUR WEBSITE

Copyright © 2023 Innomed-Europe LLC, All rights reserved.

Sie erhalten diesen Newsletter da Sie als Kunde bei Innomed-Europe registriert sind.

Our mailing address is:
Innomed-Europe LLC
Alte Steinhauserstrasse 19
6330 Cham Switzerland

Want to change how you receive these emails?
You can [update your preferences](#) or [unsubscribe from this list](#).