



Innomed eNews – Juni 2024

Angled Hip Capsule



Product No: 1767
Overall Length: 8.2" (20,8 cm)
Jaw Angle: 17°

Chandran Cannulated Dilator / Sizer for Reconstruction

Designed for dilating and sizing the bony tunnel during ACL reconstruction



Designed by Rama E. Chandran, MD


Can also be used for sizing the tenodesis screws.



Product No: 8204
Overall Length: 9" (22,9 cm)
Handle Width: 3.5" (8,9 cm)
Cannulated for Wires Up To: .062" (1,6 mm)
Sizer Markings: From 4 to 9 mm

DDA Posterior Retractor

A posterior retractor designed with a square tip and larger curvature



Product No: 6571
Overall Length: 9.75" (24,8 cm)
Blade Width: .9" (2,3 cm)
Blade Depth: 4.35" (11 cm)

Direct Anterior Angled Curette

Angled design to better access the medullary canal during anterior total hip surgery



Product No: C1031
Overall Length: 12.5" (31,8 cm)
Handle Length: 6" (15,2 cm)
Shaft Bend Angle: 30°
Curette Cup Dimensions: 8 x 14 mm

Fracture Reduction Clamp Double Curve with Speedlock

Designed with a speedlock clamp and shouldered tips to reduce bone fractures



Product No: 1755-01
Overall Length: 7.9" (20 cm)



Mini Lexer Gouges

Can be used to help remove bone from around screw heads or broken screws



Product No's:

2022-02 [Mini Lexer Gouge – 4 mm]
Overall Length: 7" (17,8 cm)
Gouge Width: 4 mm

2022-03 [Mini Lexer Gouge – 6 mm]
Overall Length: 7" (17,8 cm)
Gouge Width: 6 mm

2022-04 [Mini Lexer Gouge – 10 mm]
Overall Length: 7" (17,8 cm)
Gouge Width: 10 mm


Modified Humeral Head Reactors

Designed to help lever and displace the proximal humerus posteriorly




Product No's:

2436 [Deep]
Overall Length: 10" (25,4 cm)
Blade Width: 26 mm



2437 [Shallow]
Overall Length: 11" (30 cm)
Blade Width: 26 mm



Für weitere Fragen zu unseren Produkten stehen wir Ihnen gerne zur Verfügung!
Do not hesitate to contact us for further information about our products!

Tel. (for EU states): +49 (0) 7720 46110 60

Tel. (outside the EU): +41 (0) 41 740 67 74

For orders/offers: orders@innomed-europe.com

For loaners: loaners@innomed-europe.com

General Information: info@innomed-europe.com

VISIT OUR WEBSITE

Copyright © 2023 Innomed-Europe LLC, All rights reserved.

Sie erhalten diesen Newsletter da Sie als Kunde bei Innomed-Europe registriert sind.

Our mailing address is:
Innomed-Europe LLC
Alte Steinhauserstrasse 19
6330 Cham Switzerland

Want to change how you receive these emails?
You can [update your preferences](#) or [unsubscribe from this list](#).