





Innomed-INFO – August 2024

Levy Anterior Glenoid Retractors

Designed to help alleviate tension on anterior glenoid structures and the handle is designed to optionally be clamped to the drape



Designed by Jonathan Levy, MD



Product No's:

4536 [Small]
Overall Length: 10.5" (26,7 cm)
Depth from Bend: 5.875" (14,9 cm)
Blade Width: .75" (1,9 cm)
Tooth Gap: .325" (8,2 mm)

4536-01 [Medium] **NEW!**
Overall Length: 10.5" (26,7 cm)
Depth from Bend: 5.875" (14,9 cm)
Blade Width: 1" (2,54 cm)
Tooth Gap: .42" (10,6 mm)

4536-02 [Large] **NEW!**
Overall Length: 10.5" (26,7 cm)
Depth from Bend: 5.875" (14,9 cm)
Blade Width: 1.25" (3,2 cm)
Tooth Gap: .55" (14 mm)

Lonner Swan Lateral Knee Retractor

Ergonomically designed for more effective retraction when using a robotic arm, allowing for clearer views of the surgical site

Designed by Jess Lonner, MD and Martin Hyneman


The retractor can effectively protect the lateral soft tissues and the patella when resecting the tibial plateau and lateral femoral condyle.



Product No: 6651
Overall Length: 12.85" (32,6 cm)
Blade Width at Body: 1" (2,54 cm)
Blade Width Narrows To: .6" (1,5 cm)
Blade Width at Tip: .13" (3,3 mm)

Ortho Mallets

Designed to exert force on osteotomes, chisels, gouges and other instruments



Product No's:

7811 [Standard with Flat Sides]
Overall Length: 7.75" (19,7 cm)
Handle Length: 4.812 (12,2 cm)
Head Diameter: 1.5" (3,8 cm)
Head Weight: 2.25 lb. (1.02 kg)

7812 [Standard]
Overall Length: 7.75" (19,7 cm)
Handle Length: 4.812 (12,2 cm)
Head Diameter: 1.5" (3,8 cm)
Head Weight: 2.5 lb. (1.13 kg)

Bestellinformation:

Kontaktieren Sie uns unter +41 (0) 41 740 67 74

Für weitere Fragen zu unseren Produkten stehen wir Ihnen gerne zur Verfügung!

Do not hesitate to contact us for further information about our products!

Tel. (for EU states): +49 (0) 7720 46110 60

Tel. (outside the EU): +41 (0) 41 740 67 74

For orders/offers: orders@innomed-europe.com

For loaners: loaners@innomed-europe.com

General Information: info@innomed-europe.com

VISIT OUR WEBSITE

Copyright © 2023 Innomed-Europe LLC, All rights reserved.

Sie erhalten diesen Newsletter da Sie als Kunde bei Innomed-Europe registriert sind.

Our mailing address is:
Innomed-Europe LLC
Alte Steinhauserstrasse 19
6330 Cham Switzerland

Want to change how you receive these emails?
You can [update your preferences](#) or [unsubscribe from this list](#).